



The ibidi product family comprises a variety of different shapes of μ-Slides and μ-Dishes which all have been designed for high-end microscopic analysis of fixed or living cells. The high optical quality of the material is similar to that of glass, so you can perform all kinds of fluorescence experiments with uncompromised resolution and choice of wavelength. The μ-Slide 8 well is an array of 8 square fields where cells can be cultivated and subsequently investigated with microscopical methods. It is intended for optimization of experimental parameters like antibody dilution, seeding density or most effective drug concentration.

## Material

ibidi μ-Slides and μ-Dishes consist of a plastic with highest optical quality. The material exhibits extremely low birefringence and autofluorescence, both similar to that of glass. It is not possible to detach the bottom from the upper part. The μ-Slides and μ-Dishes are not autoclavable since they are temperature stable up to 80°C/175°F only. Please note that gas exchange between the medium and incubator's atmosphere occurs partially through the plastic bottom which should not be covered.

### Optical properties ibidi standard bottom

|                                 |                    |
|---------------------------------|--------------------|
| Refractive index $n_D$ (589 nm) | 1.52               |
| Abbe number                     | 56                 |
| Thickness                       | No. 1.5 (180 μm)   |
| Material                        | microscopy plastic |

## Geometry

The μ-Slide 8 well provides standard slide format according to ISO 8037/1.

### Geometry of μ-Slide 8 well

|                                       |                     |
|---------------------------------------|---------------------|
| Number of wells                       | 8                   |
| Dimensions of wells (w × l × h) in mm | 9.4 × 10.7 × 6.8    |
| Growth area per well                  | 1.0 cm <sup>2</sup> |
| Recommended filling volume per well   | 300 μl              |
| Total height with lid                 | 8 mm                |
| Bottom matches coverslip              | No. 1.5             |

## μ-Slide surfaces

Depending on your cells and special application you will need μ-Slides with different surfaces. If you do not need any special adhesion molecules for your application the best choice will be ibiTreat, a tissue culture treated surface.

<sup>1</sup>Collagen IV, BD Cat.-Nr. 35 6233, Fibronectin, BD Cat.-Nr. 354008, Poly-L-Lysin, Sigma Cat.-Nr. P4832, Poly-D-Lysin, BD Cat.-Nr. 354210

We provide precoated μ-Slides with adhesion substrates like Collagen IV, Fibronectin, Poly-L-Lysin, and Poly-D-Lysin. Such adhesion substrates have been shown to stimulate adhesion and growth of various cell lines in μ-Slides. Only high quality substrates are used <sup>1</sup>.

The uncoated μ-Slide is manufactured from hydrophobic plastic. For cultivation of most cell lines it is indispensable to treat the uncoated μ-Slide with biopolymers which mediate cell adhesion and growth.

## Coating your μ-Slide 8 well

The uncoated μ-Slide must be coated to promote cell adhesion. If you like to establish a certain coating for your demands we recommend testing your coating procedure on uncoated and ibiTreat μ-Slides, since we have observed that some biomolecules adhere differently to hydrophobic or hydrophilic plastic surfaces.

- Prepare your coating solution according to the manufacturer's specifications or reference.
- Apply 300 μl per well and leave at room temperature for at least 30 minutes.
- Aspirate the solution and wash with ultra-pure water. Let dry at room temperature.

## Seeding cells

- Trypsinize and count cells as usual. Dilute the cell suspension to the desired concentration. Depending on your cell type, application of a 4–9 × 10<sup>4</sup> cells/ml suspension should result in a confluent layer within 2–3 days.
- Apply 300 μl cell suspension into each well of the μ-Slide. Avoid shaking as this will result in inhomogeneous distribution of the cells.
- Cover reservoirs with the supplied lid. Incubate at 37°C and 5% CO<sub>2</sub> as usual.

## Instructions

### μ-Slide 8 well

Undemanding cells can be left in their seeding medium for up to three days and grow to confluence there. However, best results might be achieved when the medium is changed every 1–2 days. Carefully aspirate the old medium and replace it by 300 μl/well fresh medium.

#### Tip:

As you may know from 96 well plates, the bent meniscus at the air–liquid interphase in small open wells destroys the phase contrast effect of your microscope image. To avoid this problem, we recommend using our channel Slides such as the μ-Slides I and VI or the open μ-Slide 18 well flat.

### Preparation for cell microscopy

To analyze your cells no special preparations are necessary. Cells can be observed live or fixed directly in the μ-Slide on an inverted microscope. You can use any fixative of your choice. The μ-Slide material is compatible with a

variety of chemicals, e.g. Acetone or Methanol. Further specifications can be found at [www.ibidi.com](http://www.ibidi.com). Due to the thin bottom of only 180 μm, high resolution microscopy is possible.

### Immersion oil

When using oil immersion objectives, only the immersion oils specified in the table may be used. The use of different oil can lead to damages of the plastic material and the objective.

| Company  | Product                         | Ordering number   |
|----------|---------------------------------|-------------------|
| Cargille | type DF, Formula Code: 1261     | (Cargille) 16242  |
| Zeiss    | 518 F                           | (Zeiss) 444960    |
| Olympus  | 50CC                            | (Olympus) 35506   |
| Nikon    | 50 CCM DF                       | (Nikon) MXA 20351 |
| Leica    | immersion oil, low fluorescence | (Leica) 11513859  |

### μ-Slide 8 well family

The μ-Slide 8 well family is available with different surfaces. See table below for choosing your μ-Slide 8 well.



| Ordering number | Treatment or Coating    | characteristics                     |
|-----------------|-------------------------|-------------------------------------|
| 80826           | ibiTreat, sterile       | hydrophilic, tissue culture treated |
| 80822           | Collagen IV, sterile    | protein coating                     |
| 80823           | Fibronectin, sterile*   | protein coating                     |
| 80824           | Poly-L-Lysine, sterile  | biopolymer coating                  |
| 80825           | Poly-D-Lysine, sterile* | biopolymer coating                  |
| 80821           | uncoated, sterile       | hydrophobic                         |

\* available on request only

**Selected cell tests on different surfaces**

Many eukaryotic and bacterial cells have been tested by ibidi on the different surfaces of the µ-Slides. A variety of other cell lines like COS, CHO, HepG2, and NIH 3T3 were successfully grown by our customers.

|                                 | ibiTreat  | Collagen IV | Fibronectin | Poly-L-Lysin   | Poly-D-Lysin | uncoated       |
|---------------------------------|-----------|-------------|-------------|----------------|--------------|----------------|
| HUVEC                           | excellent | good        | excellent   | no cell growth | not done     | no cell growth |
| Rat1                            | excellent | excellent   | excellent   | excellent      | excellent    | poor           |
| HT1080                          | excellent | excellent   | excellent   | excellent      | not done     | poor           |
| HeLa                            | excellent | excellent   | excellent   | excellent      | not done     | poor           |
| Neuro2A                         | excellent | excellent   | excellent   | excellent      | excellent    | poor           |
| PC12                            | good      | excellent   | excellent   | excellent      | excellent    | no cell growth |
| <i>Dictyostelium discoideum</i> | not done  | excellent   | not done    | not done       | not done     | excellent      |
| <i>Escherichia coli</i>         | excellent | not done    | not done    | excellent      | not done     | excellent      |

HUVEC = Human Umbilical Vein Endothelial Cells

Rat1 = Rat Fibroblast

HT1080 = Human Fibrosarcoma

HeLa = Human Cervix Adenocarcinoma

Neuro2A = Mouse Neuroblastoma

PC12 = Rat Pheochromocytom

*Dictyostelium discoideum* = strain wild type AX-2

*Escherichia coli* = strain MDG131

**For research use only!**

Further technical specifications can be found at [www.ibidi.com](http://www.ibidi.com). For questions and suggestions please contact us by e-mail [info@ibidi.de](mailto:info@ibidi.de) or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany.

© ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.